

Tyler Chang<sup>1</sup>, Jakob Elias<sup>1</sup>, Stefan Wild<sup>2</sup>, Santanu Chaudhuri<sup>1,3</sup> and Joseph Libera<sup>1</sup>


<sup>1</sup>Argonne National Laboratory <sup>2</sup>Lawrence Berkeley National Laboratory <sup>3</sup>University of Illinois Chicago

## Our (Big) Goals

- Design a software framework for *self-driving labs*
- Accelerate discovery via intelligent experimentation
- Democratize lab-work by building open-source tools


## Streaming data from multiple sources

- How to collect and analyze data?
- MDML is a platform for streaming, analyzing, and logging experiment and simulation data



## Our framework and software stack

- MDML gives us access to *heterogeneous data from laboratory sources*
- ParMOO gives us modular/customizable *modeling, embedding, and solvers*




Build & deploy custom solvers for computational and experimental problem!

## Experiment Results

- **Want to maximize TFMC production at high temperatures**
- High temperatures trigger a side-reaction and produces byproduct (TFE)
- 15-pt Latin hypercube, Gaussian RBF surrogate, L-BFGS-B optimizer
- 3 scalarizations per batch, sorted by temp
- evaluated batch on CFR, TFMC and TFE peaks recorded by NMR

Results after 41 experiments steered by our solver




Get this code!

```
git clone https://github.com/parmoo/cfr-materials
pip install REQUIREMENTS.txt
```

## Next Steps

Need to handle more complex design spaces:

- Generative AI for embeddings
- Trust-region descent methods
- Subspace iterations
- Custom surrogate models
- Structure-exploiting optimizers




(Top) using custom embedders to optimize in latent space

(Bottom) exploiting problem structure using composite objectives



## Model-based optimization

- ParMOO (multiobjective optimization) library is used to implement an active learning loop



## Example: TFMC Manufacturing Conditions

Optimize the production of TFMC via a known reaction...



... in LabVIEW automated CFR and measured via NMR

Left to right:  
PC running LabVIEW CFR and feed NMR spectroscope

Design variables and bound constraints for experiment

Parameter	Lower bound	Upper bound
Temperature (degrees C)	40	150
Reaction time (seconds)	60	300
Equivalence ratio (no units)	0.9	2

This material was based upon work supported by the U.S. Department of Energy, Office of Science, Office of Advanced Scientific Computing Research, Applied Mathematics and SciDAC programs under Contract Nos. DE-AC02-05CH11231 and DE-AC02-06CH11357 and by the Argonne LDRD program.